

# TILDA

## - Tool for Interactive Learning and Daily Assistance

TILDA is a web based learning tool, which through interactive training courses makes it possible to secure quality and quality control education within the healthcare system. TILDA is developed in collaboration with the Intensive Care unit/Anaesthesia unit at the University Hospital MAS (Malmö).

### How does TILDA work?

TILDA consists of two main units. With the administrative unit of TILDA its possible to create interactive training of medical technology, methods for operator certification and specific training on your unit/department. Training courses are fast and easy to create by means of TILDA's training course templates. Users taking part in the interactive training courses are guided through theoretical, practical parts and questions in a educational manner. The theoretical, practical parts and the questions are finished with a test and can be sandwiched with practical reviews.

TILDA is accessible night and day and each user gets a personally adapted account, this makes it possible for the user, to work with the tool, on its own, whenever possible.


### Quality control

TILDA is quality controlling the training by personal certificates. The personal certificate is a conformation that the user has taken part in the training course and utilised sufficient knowledge. In this way, the personal certificate, creates quality control and safety for the caregivers and the caretakers.

The personal certificates makes it possible for administrative users to get a general view of the results of the training course of the unit/department. The general view can be presented as clear statistics and may be read directly on the screen or printed which makes it possible to follow up the activities

### Get a quick start

Because of TILDA's clear and useful layout, may experienced as well as inexperienced users use the tool in a fast and effective way. TILDA is flexible and can be tailor-made to fit the needs and desires of a unit/department. TILDA can quickly be implemented to a large number of users.



### Contact us

Contact us for more information about TILDA at [info@interactit.se](mailto:info@interactit.se). For more information about how TILDA is used at an unit/department please contact Ms. Sue Harden, Spec Nurse at the Intensive Care Unit at the University hospital MAS (Malmö), [sue.harden@skane.se](mailto:sue.harden@skane.se)

#### 1. Register and login

At the registration a user specifies his or hers personal information, TILDA automatically adapts to the user's occupational status.

#### 2. Choose training course

The user chooses an interactive training course, which can comprise everything, from a device to a working method.



#### 3. Training course

The user receives the training course that is displayed, e.g. as text information, movie sequences, slide shows and/or instruction manuals. The training course is completed by interactive questions.

#### 4. Certificate of knowledge

TILDA indicates that a training course is completed and the user receives a certificate of a successful training course.